

# Dissolution Method Development for Assessing Amorphous Solid Dispersions

Manjeet Pimparade, Lawrence Martin, Manish Rane and Ali Rajabi-Siahboomi

Colorcon, Inc., Harleysville, PA, USA, [www.colorcon.com](http://www.colorcon.com)

AAPS  
Poster Reprint 2018

## Purpose

The purpose of this study was to develop a discriminatory dissolution method to assess different formulations of itraconazole amorphous solid dispersions (ASD). Various dissolution conditions such as; media composition, volume and hydrodynamics were investigated.

## Methods

ASDs of itraconazole (ITR) using hypromellose acetate succinate (AFFINISOL™ HPMCAS 716G, International Flavors and Fragrances Inc., USA) were produced by hot melt extrusion in 1:1, 1:2 and 1:3 ratios using a twin-screw extruder (Pharma 11, Thermo Fisher)<sup>1</sup>. The extrusion was carried out at 3-7 g/min feed rate, 100 rpm screw speed and target process temperature of 170°C. The extruded strands were air-cooled, pelletized and milled into a powder, then stored in double whirl packs at room temperature for evaluation within one week of processing.

Dissolution testing was carried out on ITR+HPMCAS 716G (1:3 ratio) ASD samples, in USP Type II (paddle) apparatus at 37°C, under variable conditions such as media composition (pH 1.2 to pH 6.8), presence of sodium lauryl sulphate (SLS) surfactant at concentrations of 0 to 2.0%w/w, media volume (500 and 1000 mL) and stirring speed (50, 75, 100 rpm) (Table 1). Formulations containing different drug to polymer ratio were subjected to selected dissolution methods to check the suitability to provide adequate discrimination.

**Table 1: Dissolution Conditions Studied During Development of Discriminatory Dissolution Test for ASD**


| Condition                  | Media pH            | Surfactant level (%) | Volume (mL) | Paddle Speed (rpm) | Dose of API (mg) |
|----------------------------|---------------------|----------------------|-------------|--------------------|------------------|
| Effect of pH               | pH 1.2 <sup>^</sup> | ----                 | 500         | 75                 | 25               |
|                            | pH 4.5 <sup>#</sup> | ----                 | 500         | 75                 | 25               |
|                            | pH 5.5 <sup>*</sup> | ----                 | 500         | 75                 | 25               |
|                            | pH 6.0 <sup>*</sup> | ----                 | 500         | 75                 | 25               |
|                            | pH 6.4 <sup>*</sup> | ----                 | 500         | 75                 | 25               |
| Effect of surfactant (SLS) | pH 6.8 <sup>*</sup> | 0.10                 | 500         | 75                 | 25               |
|                            | pH 6.8 <sup>*</sup> | 0.25                 | 500         | 75                 | 25               |
|                            | pH 6.8 <sup>*</sup> | 0.50                 | 500         | 75                 | 25               |
|                            | pH 6.8 <sup>*</sup> | 1.00                 | 500         | 75                 | 25               |
|                            | pH 6.8 <sup>*</sup> | 1.50                 | 500         | 75                 | 25               |
| Effect of media volume     | pH 6.8 <sup>*</sup> | ----                 | 500         | 75                 | 100              |
|                            | pH 6.8 <sup>*</sup> | 1.00                 | 500         | 75                 | 100              |
|                            | pH 6.8 <sup>*</sup> | ----                 | 1000        | 75                 | 100              |
| Effect of paddle speed     | pH 6.8 <sup>*</sup> | ----                 | 1000        | 50                 | 100              |
|                            | pH 6.8 <sup>*</sup> | ----                 | 1000        | 75                 | 100              |
|                            | pH 6.8 <sup>*</sup> | ----                 | 1000        | 100                | 100              |

## Results

### Effect of Media pH

Figure 1 shows the effect of media pH on drug release profiles using similar media volume (500 mL) and paddle stirring speed (75 rpm). Drug release (about 80%) only occurred in dissolution media with pH 6.4 and above. This behavior may be attributed to the enteric nature of HPMCAS polymer. Based on this data, all subsequent dissolution studies were conducted in phosphate buffer pH 6.8.


**Figure 1: Effect of Media pH on ITR Release (25 mg dose) from ITR:HPMCAS ASD, 1:3 Ratio**



**Effect of Surfactant (SLS) Concentration**

ITR + HPMCAS 716G dispersion was subjected to dissolution testing in sodium phosphate buffer pH 6.8 containing SLS, as surfactant, at various concentration (0.1 to 2.0%). Drug release in 0.1% SLS level was high initially, and then started to drop. All other concentrations showed an increase in drug release, up to 1 hour (Figure 2). None of the SLS concentrations achieved  $\geq 80\%$  drug release till 1 hour, with exception of media containing 1.5% SLS. Drug release in 1 to 2% SLS were similar at the end of 1.5 hour time point (data not shown).


**Figure 2: Effect of Surfactant (SLS) Concentration on ITR Release (25 mg dose) from ITR:HPMCAS ASD, 1:3 ratio**



**Effect of Media Volume**

Figure 3 shows that drug release from itraconazole ASD (100 mg dose) in low (500 mL) or high volume (1000 mL) phosphate buffer pH 6.8 was similar. However, the presence of 1% SLS gave slightly different results in two media volumes.


**Figure 3: Effect of Media Volume on ITR Release from ASD (100 mg dose) in Phosphate Buffer pH 6.8**



### Effect of paddle speed

Initial dissolution testing was carried out at 75 rpm, based on findings from the literature review.<sup>2</sup> Further study was conducted at low (50 rpm) and high (100 rpm) paddle speed. As seen in Figure 4, dissolution was lower at 50 rpm, whereas, dissolution was comparable at 75 and 100 rpm paddle speeds. Higher paddle speed allowed ASD to be uniformly dispersed in the dissolution media within 2 minutes.

**Figure 4: Effect of Paddle Speed on ITR Release from ASD**




### Determination of suitable dissolution method to provide discrimination between formulations

Based on the earlier results, formulations of ITR : HPMCAS solid dispersions at various drug : polymer ratio, were subjected to 1000 mL phosphate buffer pH 6.8, with or without 1% SLS, USP apparatus II (paddle) at 75 rpm and 37°C. Release profiles of ITR from the various ASD formulations showed adequate discrimination when subjected to phosphate pH 6.8 dissolution, as shown in Figure 5a. However, incorporation of 1% SLS in dissolution media failed to provide discrimination between various formulations. There was a definite improvement in the dissolution of

crystalline ITR, or its physical mixture with HPMCAS, in media containing 1% SLS. As reported previously, surfactants may interact with certain polymers during solubilization of API from amorphous solid dispersions, which may result in altered drug solubility, drug dissolution and drug absorption.<sup>3-4</sup>

**Figure 5: ITR Release from ASDs with Different Drug:Polymer Ratio in (a) Phosphate Buffer pH 6.8 and (b) Phosphate Buffer pH 6.8 with 1% SLS.**



## Conclusions

Media composition was shown to have a major effect on dissolution of ITR from ASD containing HPMCAS. The dissolution method that provided discrimination amongst different formulations is 1000 mL phosphate buffer pH 6.8 in USP Type II (paddle) at 75 rpm. This method was selected for further formulation evaluation. Discriminating robust dissolution methods are critical during formulation screening and drug development.

## References

1. Martin, L., Pimparade, M. et al. Fundamental evaluation and characterization of itraconazole solid dispersions prepared by hot melt extrusion. CRS 2018.
2. Pack, B., Babayan, Y. et al. Development of an in vivo-relevant drug product performance method for an amorphous solid dispersion. J Pharm Biomed Analysis, 2017, 142, 307-314.
3. Deshpande, T., Shi, H. et al. Investigation of polymer-surfactant interactions and their impact on itraconazole solubility and precipitation kinetics for developing spray dried amorphous solid dispersions. Mol. Pharmaceutics, 2018, 15(3), 962-974.
4. Chen, Y., Wang, S. et al. Sodium lauryl sulfate competitively interacts with HPMC-AS and consequently reduces oral bioavailability of posaconazole / HPMC-AS amorphous solid dispersion. Mol. Pharmaceutics, 2016, 13(8), 2787-2795

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade, copyright, patent or other rights held by any third person or entity when used in the customer's application.

For more information, contact your Colorcon representative or call:

|                                  |  |                                  |                          |                          |
|----------------------------------|--|----------------------------------|--------------------------|--------------------------|
| North America<br>+1-215-699-7733 | Europe/Middle East/Africa<br>+44-(0)-1322-293000 | Latin America<br>+54-1-5556-7700 | India<br>+91-832-6727373 | China<br>+86-21-61982300 |
|----------------------------------|--|----------------------------------|--------------------------|--------------------------|

You can also visit our website at [www.colorcon.com](http://www.colorcon.com)



© BPSI Holdings LLC, 2018.

The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

AFFINISOL™ is a trademark of International Flavors and Fragrances Inc. or its affiliates. © 2021 IFF. All rights reserved

All trademarks, except where noted, are property of BPSI Holdings, LLC.