

Evaluation of Film Coating Weight Uniformity, Tablet Progression and Tablet Transit Times in a High Throughput Continuous Coating Process

Chris Neely, Charles Cunningham and Ali Rajabi-Siahboomi
www.colorcon.com

The developmental film coating from this study was commercially launched in April 2016 as:

Opadry® QX
Quick and FlexIBLE Film Coating System

AAPS
Poster Reprint 2015

Purpose

Compared to batch coating processes, continuous film coating processes allow for fast coating times and have been shown to provide excellent color uniformity due to a shallower tablet bed depth and increased frequency of tablet presentations to the spray zone.¹ It has also been shown that high solids concentration ($\geq 20\%$) coating systems may be well suited to maximize the efficiency of continuous processes.²

Interest in continuous coating for pharmaceutical products is increasing and there is a need to more fully understand how these processes compare to traditional batch coating processes in terms of actual coating weight uniformity rather than just color uniformity. This study evaluates the effect of tablet residence time and uniformity of tablet progression on coated tablet weight variation in a high throughput continuous coating process. In addition, the coating uniformity values obtained in the continuous coating process were compared to coating uniformity in a typical production scale batch coating process.

Methods

Trials were conducted in a 24" diameter continuous tablet coater (Model CTC, Thomas Engineering, USA) equipped with 22 spray guns (ABC type, Schlick, DE) (Figure 1).

Figure 1. Exterior and Interior Views of Thomas CTC



Colorcon placebo tablets (300 mg, 10 mm, round with debossed logo) were used as the substrate. An innovative, pigmented coating formulation developed by Colorcon based on Kollicoat IR (polyvinyl alcohol-polyethylene glycol graft copolymer) (BASF, Florham Park, NJ) was used in the study. The formulation was prepared at both 20% or 25% solids concentration and applied to a target weight gain (WG) of 3.0% at tablet throughput rates of both 680 kg/hour or 850 kg/hour using the process parameters shown in Table 1.

Table 1. Coating Process Parameters

CTC Parameters	Trial 1	Trial 2
Inlet air temperature (°C)	75	72
Inlet airflow rate (cfm / m ³ /hr)	6700 / 11383	6600 / 11213
Target exhaust temperature (°C)	47	47
Target bed temperature (°C)	45	45
Solids concentration (%)	20	25
Dispersion spray rate (g/min)	1700	1700
Atomizing air pressure (psi/bar)	35 / 2.4	35 / 2.4
Pattern air pressure (psi/bar)	30 / 2.1	30 / 2.1
Pan speed (rpm)	12	12
Resultant EEF*	3.26	3.25
Tablet feed rate (kg/hr)	680	850
Target resident pan load (kg) <small>(quantity of tablets in the pan at any given time)</small>	120	120
Coating weight gain (%)	3.0	3.0

* Environmental Equivalency Factor (Thomas Engineering TAAC Program)

Individually marked tracer tablets were used to determine tablet transit times and coating weight variation. Each tracer tablet was marked with a unique letter and number. The marked tablets were dried to constant weight in a 50°C oven. After drying, the weight of each individual tablet was recorded. At the start of each trial, groups of 100 marked tablets were added to the in-feed side of coating process every 30 seconds for 2.5 minutes for a total of 500 marked tablets in each trial. The time of addition for each group of marked tablets was recorded.

As the coating progressed, all coated tablets were captured from the discharge side of the pan in individual bags at 30 second intervals. After completion of each coating trial, the tracer tablets were sorted from each bag and the time of discharge for each tablet was recorded. This allowed tablet transit times through the coater to be determined with an accuracy of ± 30 seconds. The coating formulation used a low level of Yellow #10 pigment to allow for visual assessment of color uniformity but still allow for the markings on the tracer tablets to be clearly seen (Figure 2).

Figure 2. Example of Marked Tracer Tablets Before and After Coating



The recovered tracer tablets were then re-dried to constant weight and re-weighed to determine the actual coating WG applied to each tablet as well as coating weight variation amongst each set of marked tablets. Drying the tablets to constant weight before and after coating ensured that any moisture gain or loss as a result of the coating process would not impact the accuracy of coating weight gain determination.

Comparative trials were conducted using marked tablets in a 48" fully perforated batch pan at the same 3.0% target WG, coating solids concentrations, and appropriate process scale-out factors to maintain similar thermodynamic conditions (EEF) and linear pan velocity. The process conditions used for these trials are shown in Table 2.

Table 2. Process Conditions Used for 48" Batch Coating Trial

48" Batch Pan Parameters	Trial 1	Trial 2
Inlet air temperature (°C)	62	62
Inlet airflow (cfm / m ³ /hr)	1800 / 3058	1800 / 3058
Target exhaust temperature (°C)	48	47
Target bed temperature (°C)	47	45
Solids concentration (%)	20	25
Dispersion spray rate (g/min)	400	400
Atomizing air pressure (psi/bar)	25 / 1.7	25 / 1.7
Pattern air pressure (psi/bar)	25 / 1.7	25 / 1.7
Pan speed (rpm)	6	6
Resultant EEF	3.23	3.22
Pan load (kg)	130	130
Coating weight gain (%)	3.0	3.0

Results

Coated Tablet Appearance – Continuous Coating Process

All coated tablets were uniform in appearance with no visually apparent defects. Excellent logo definition was noted for all samples (Figure 3).

Figure 3. Coated Tablet Appearance



Tablet Residence and Transit Times

Tablet transit times through the continuous coater averaged 15 min at 680 kg/hr and 14 min at 850 kg/hr respectively. For each trial 50% of all marked tablets exited the coater within a span of ~2-3 minutes while the remaining 50% of the marked tablets exited the coater in an elapsed coating time as short as 11.5 minutes.

In the case of a few tablets, 27 minutes of coating time elapsed. Statistical representations of tablet transit times are shown in Figures 4 and 5.

Figure 4. Tablet Progression Through the CTC at 680 kg/hr (20% Solids Concentration)






Figure 5. Tablet Progression Through the CTC at 850 kg/hr (25% Solids Concentration)



Coated Tablet Weight Gain and Variability – Continuous Coating Process

For Trial 1, coated tablet WG averaged 3.0% as targeted. For Trial 2, the average WG was slightly below target at 2.7% WG. It is suspected that there was a slight increase in the resident pan load at the higher throughput rate of 850 kg/hr compared to the previous trial at 680 kg/hr. Surprisingly, it was found that coated tablet weight gain was independent of residence time despite some tablets remaining within the coater for a substantially longer period of time than others. It is likely that some tablets appear more frequently in front of the spray zone while others appear less frequently, similar to typical tablet movement in traditional batch coating. This relationship between elapsed coating time and coating weight gain is shown with 4 sets of marked tablets in Figure 6.

Figure 6. Effect of Tablet Residence Time on Coating Weight Gain in the CTC



Coated tablet weight variation ranged from 15.25% relative standard deviation (RSD) to 22.37% RSD. Increasing throughput rates and higher coating solids concentration resulted in increased weight variability as shown in Table 3.

Table 3. Coating Weight Gain and Variability

Trial	Tablet Set	Throughput Rate (kg/hr)	Solids Conc. %	Avg. Weight Increase (mg) (n=100)	St. Dev. (mg)	Actual Weight Gain (%)	% RSD
1	A	680	20	10.57	1.74	3.02	16.49
	B			10.75	1.60	3.07	14.93
	C			10.70	1.62	3.06	15.25
	D			10.39	1.72	2.97	16.66
	E			10.46	1.67	2.99	16.14
2	F	850	25	9.49	1.82	2.71	19.11
	G			9.21	2.05	2.62	22.19
	H			9.49	1.99	2.71	21.22
	J			9.28	1.95	2.65	20.82
	K			9.15	2.02	2.61	22.37

Coating Uniformity Comparison with Traditional Batch Coater

Reference trials in the 48" batch pan exhibited higher coated tablet weight variation than trials conducted in the CTC. This is especially significant given the much longer actual coating time in the batch coater vs. the CTC. A comparison between coating uniformity in the CTC vs. the traditional batch coater is shown in Table 4.

Table 4. Coating Time and Coating Uniformity Comparison – CTC vs. Batch Coating

Parameters	CTC	48" Batch	CTC	48" Batch
Solids conc. (%)	20		25	
Pan load (kg)	120 (resident)	130	120 (resident)	130
Total coating time (min)	15 (average residence time)	49	14 (average residence time)	39
Coated tablet weight variability (% RSD)	15.90	18.74	21.10	23.90

Conclusions

Coating variability and tablet movement in a production scale continuous coating process were quantified. Tablets progressed through the continuous process with relatively small variability in transit times and exhibited superior uniformity compared to tablets coated in a traditional batch coater. Higher productivity coating systems ($\geq 20\%$ solids) are better suited for continuous processes where coating weight variability is reduced relative to traditional batch coating processes. Further work is ongoing to assess the effects of changing pan speeds and even higher solids concentration coatings in continuous coating processes.

Acknowledgements

The authors are grateful to Thomas Engineering Inc. for the generous use of their equipment, assistance with trials and coating expertise. We would also like to thank Joseph Quinn of Colorcon Inc. for his careful preparation and analysis of the marked tablets used in these trials.

References

- Cunningham C, Hansell J, Nuneviller III F, Rajabi Siahboomi AR (2010). Evaluation of recent advances in continuous film coating processes. Drug Development and Industrial Pharmacy (36) 2 :227-233.
- Neely C, Cunningham C, Rajabi Siahboomi AR (2014). Evaluation of a high productivity developmental aqueous coating formulation in a high throughput continuous coating process. Poster presentation, AAPS Annual Meeting and Exposition, San Diego, CA.

This submission is an Encore presentation originally presented on March 16-17, 2015 at the 50th Annual AAPS Arden Conference on Continuous Processing in Baltimore, MD.

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

For more information, contact your Colorcon representative or call:

North America
+1-215-699-7733

Europe/Middle East/Africa
+44-(0)-1322-293000

Asia Pacific
+65-6438-0318

Latin America
+54-1-5556-7700

You can also visit our website at www.colorcon.com



©BPSI Holdings, LLC 2015

All trademarks, except where noted, are property of BPSI Holdings LLC. The information contained in this document is proprietary to Colorcon, Inc. and may not be used or disseminated inappropriately.